

Blue to the sky. Orange to the thigh.

How to use EpiPen® and EpiPen® Jr (epinephrine) Auto-injectors.

Remove the EpiPen® Auto-Injector from the carrier tube and follow these 2 simple steps:



- Hold firmly with orange tip pointing downward.
- Remove blue safety cap by pulling straight up. Do not bend or twist.



- Swing and push orange tip firmly into mid-outer thigh until you hear a 'click'.
- Hold on thigh for several seconds.



Built-in needle protection

- After injection, the orange cover automatically extends to ensure the needle is never exposed.



After administration, patients should seek medical attention immediately or go to the emergency room. For the next 48 hours, patients must stay within close proximity to a healthcare facility or where they can call 911.

For more information visit EpiPen.ca.

Indications and Clinical Use:

- EpiPen® and EpiPen® Jr Auto-Injectors are indicated for the emergency treatment of anaphylactic reactions in patients who are determined to be at increased risk for anaphylaxis, including individuals with a history of anaphylactic reactions.
- Select dose according to patient's body weight.
- Inject into the outer thigh.
- Following treatment of anaphylaxis, the patient must call 911.

There are no absolute contraindications to the use of epinephrine in a life-threatening allergic situation.

For More Information:

Please consult the prescribing information at http://www.pfizer.ca/en/our_products/products/monograph/312 for more information relating to adverse reactions, drug interactions, and dosing information, which has not been discussed in this piece. The prescribing information is also available by calling 1-800-463-6001.

Relevant Warnings and Precautions:

- Avoid use of epinephrine in patients with cardiogenic, traumatic, or hemorrhagic shock; cardiac dilation; cerebral arteriosclerosis; organic brain damage, and narrow-angle glaucoma
- Use with caution in patients with cardiac arrhythmias; coronary artery or organic heart disease; hypertension and hyperthyroidism
- Sulfite sensitivity
- Diabetes
- Parkinson's disease



© 2014 Pfizer Canada Inc., Kirkland, Quebec H9J 2M5 • Toll free: 1-877-EPIPEN1 (1-877-374-7361)
EpiPen®, EpiPen® Jr are registered trademarks of Mylan Inc. licensed exclusively by its wholly-owned affiliate,
Mylan Specialty, L.P., sub-licensee, Pfizer Canada Inc., Kirkland, Quebec H9J 2M5
D0000487988



Trusted for over 25 years.